

## APN-1 Apnoea Monitor

Utilising advanced digital techniques the APN-1 Apnoea monitor is designed to monitor breathing patterns but ignore unwanted movements such as the heartbeat.

Each breath is indicated by the APN-1 Apnoea monitor with a flashing green light.

Should breathing cease for more than 20 seconds the local alarm will sound.

The control box is placed on the bedside table, with the sensor pad being placed under the mattress. For best performance the monitor should be used with a foam type mattress.



### Features

### Benefits

#### No skin contact

The APN-1 monitors breathing movements with an under-the-mattress sensor pad.

#### Breathing Indications

As each breath is detected a green light will flash and a brief 'tick' sounds.

#### Apnoea Alarm

If breathing movements stop, there will be a pre-warning bleep after 15 Seconds then the main alarm sounds at 20 seconds.

#### Battery Low Warning

The Monitor will provide visible and audible indication prior to the batteries needing replacement.

#### Signal Level Indication

The Monitor will indicate if the signal level from the sensor pad is too low for reliable detection, the user may correct any problems and avoid false alarms.

#### Breathing Rate Alarm

The Monitor will produce a different alarm sound if the Baby's breathing falls below or rises above 'normal' rates. The alarm light will also flash slowly for slow breathing and rapidly for fast breathing.

#### Remote Alarm Option

To extend the range of the monitor, a transmitter and remote receiver can be purchased, this allows the alarm to be transmitted to another room up to 30m away.