

Datasheet

BC2 Pneumatic BreathCall Switch

BreathCall Pneumatic Air Switch

- Blow activated air switch
- Blow into the mouthpiece to activate the call point
- Bedside Table clamp
- Disposable filters for infection control
- Laboratory tested for use in an oxygen enriched atmosphere
- Suitable for children and adults
- Links to directly to Intercall and most other Nursecall systems*



*Standard contact type: Normally open ¼" Mono Jack

Features

Blow Switch

The BreathCall unit can be used as an aid to communication for persons with limited dexterity that, are unable to use a standard call lead.
eg. partially paralysed individuals and those with severe arthritic and rheumatic conditions.

Easy Positioning

The flexible tube means that the BreathCall unit can be positioned to suit the individual and can be easily transferred from one position to another and fixed in position with the table clamp.

Activation

Requires only the lightest blow to activate, the BreathCall unit will then trigger the Alarm system or Nursecall system to alert the carer that assistance is required.

Operation

This is a N/O vacuum switch going closed on activation, power is not necessary. It can accept up to 24 volts AC or DC at .33amps.