

Product Datasheet

PC1

Pneumatic Presscall Switch

- Soft-touch air bulb activated
- Laboratory tested for use in an oxygen enriched atmosphere
- Equipped with stainless steel bed clip
- Suitable for children and adults
- Links to directly to Intercall and most other Nursecall systems*

*Standard contact type: Normally open ¼” Mono Jack



Features

Pneumatic Touch

Benefits

The Presscall unit can be used as an aid to communication for persons with limited dexterity that, are unable to use a standard call lead. For example: partially paralysed individuals and those with severe arthritic and rheumatic conditions.

Easy Positioning

The flexible tube means that the Presscall unit can be positioned to suit the individual and can be easily transferred from one position to another and fixed in position with the stainless steel bed clip.

Activation

Requires only the lightest touch to activate, the Presscall unit will then trigger the Alarm system or Nursecall system to alert the carer that assistance is required.

Operating Voltage

Up to 24 volts AC or DC